

Tissue Specific Antibodies in Early Sjögren's Disease Diagnosis

Includes the Following Biomarkers:

Salivary protein 1 antibodies - IgG, IgA & IgM;
Carbonic anhydrase VI antibodies - IgG, IgA & IgM,
& Parotid specific/secretory protein antibodies -
IgG, IgA, & IgM)

KSL Diagnostics Test Code:

308

Methodology:

Chemiluminescence

Units:

U/ml

Reference Range:

- <20 U/ml - Negative
- 20-25 U/ml - Borderline
- >25 U/ml- Positive

CPT Code:

83520 (x9)

Schedule/Turnaround Time:

Assay performed once weekly on Wednesday.
Report availability is within one week from the time
of specimen receipt.

Specimen Requirements:

Collect 5-10 ml of blood in a red top or serum
separator tube. Separate serum from cells ASAP or
within 2 hours of collection. Transfer 1 mL serum to
a Standard Transport Tube. Specimen need not be
refrigerated or frozen.

Requested Specimen Volume:

2 mL

Absolute Minimum Volume:

0.5 mL

Cause for Rejection:

Specimens other than serum. Grossly hemolyzed,
lipemic or icteric samples.

Sample Stability:

Sample is stable at ambient temperature during
shipment. If sample is stored prior to shipment, it
is stable refrigerated (2-8°C) up to five days and
frozen (-20°C or lower) up to one year.

Clinical Relevance:

Sjogren's syndrome (SS) is a systemic autoimmune
disease in which loss of salivary gland and
lachrymal gland function is associated with
hypergammaglobulinemia, autoantibody
production, mild kidney, and lung disease and
eventually lymphoma. SS involves dry eyes and dry
mouth without systemic features that may be either
primary or secondary to another autoimmune
disease, such as SLE in patients with SS diagnosed
at a late stage in their disease, after the salivary
glands and lachrymal glands are already destroyed.
Only symptomatic treatment can be offered for
abnormal lachrymal and salivary gland function.
The diagnosis for SS is currently at a crossroad with
the American College of Rheumatology providing
which requires characteristic autoantibodies
(SS-A) or minor salivary gland biopsy. Since lip
biopsies are not frequently performed in clinical
practice, there is increased emphasis placed on
autoantibodies in diagnosis. Novel antibodies
specific to tissues in salivary and lacrimal glands
were identified in 2012 by Shen et al., includes
salivary gland protein 1 (SP-1), carbonic anhydrase
VI (CA-VI) and parotid secretory protein (PSP).
Further studies have shown that the isotype
differentiation of the markers adds to the sensitivity
of diagnosis of SS. These autoantibodies occurred
earlier in the course of the disease than antibodies
to Ro or La. In addition, antibodies to SP-1, CA-
VI and PSP were found in patients meeting the
criteria for SS who lacked antibodies to Ro or La.
Furthermore, in patients with idiopathic xerostomia
and xerophthalmia for less than 2 years, 76% had
antibodies to SP-1 and/or CA6 while only 31%
had antibodies to Ro or La. Antibodies to different
isotypes (IgG, IgM & IgA) of SP-1, CA-VI and PSP
are useful markers for identifying patients with SS
at early stages of the disease or those that lack
antibodies to either Ro or La.

Selected References:

- Fox, R (2005). Sjogren's syndrome. Lancet;
366: 321–331. 27
- Shen, L. et al., (2012). Novel autoantibodies in
Sjogren's syndrome. Clinical Immunology;145,
251–255.